

**Perforated Plate Test Sieves with Round Apertures, Ø 200 mm, Stainless Steel** **Blau-Metall Specs & Prices (01/2025)**

Aperture (w)	Margin / Edge	Perf. Area (ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	200 x 32 mm			200 x 50 mm		
							Art-No.	€/Pcs.	Deliv. Time*	Art-No.	€/Pcs.	Deliv. Time*
0,50 mm	Perforated	200	1,09	0,5 mm	Rv	DIN 24041**	520050.32	135,00 €	7-14 days	520050.50	137,00 €	4-10 days
0,75 mm	Perforated	200	1,50	0,6 mm	Rv	Tolerances ISO 3310	520052.32	168,00 €	7-14 days	520052.50	170,00 €	4-10 days
1,00 mm	Perforated	200	2,00	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512100.32	151,00 €	7-14 days	512100.50	153,00 €	7-10 days
1,25 mm	Perforated	200	2,40	1,0 mm	Rv	ISO 3310-2	512103.32	151,00 €	7-14 days	512103.50	153,00 €	4-10 days
1,40 mm	Unperforated	185	2,60	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512104.32	156,00 €	7-14 days	512104.50	158,00 €	7-10 days
1,40 mm	Unperforated	185	2,60	0,8 mm	Rv	EN ISO 5223	732250.32	166,10 €	7-14 days	732250.50	168,10 €	7-10 days
1,50 mm	Unperforated	185	2,50	0,8 mm	Rv	Tolerances ISO 3310	520054.32	156,00 €	7-14 days	520054.50	158,00 €	7-10 days
1,55 mm	Unperforated	185	2,50	0,8 mm	Rv	Tolerances ISO 3310	520155.32	168,00 €	7-14 days	520155.50	170,00 €	7-10 days
1,60 mm	Unperforated	185	2,75	0,8 mm	Rv	ISO 3310-2	512105.32	156,00 €	7-14 days	512105.50	158,00 €	7-10 days
1,80 mm	Unperforated	185	3,10	0,8 mm	Rv	ISO 3310-2	512107.32	151,00 €	7-14 days	512107.50	153,00 €	7-10 days
1,80 mm	Unperforated	185	3,10	0,8 mm	Rv	EN ISO 5223	732251.32	166,10 €	7-14 days	732251.50	168,10 €	7-10 days
1,90 mm	Unperforated	185	3,30	0,8 mm	Rv	Tolerances ISO 3310	520183.32	168,00 €	7-14 days	520183.50	170,00 €	7-10 days
1,98 mm	Unperforated	185	3,97	0,8 mm	Rv	U.S. CFR	622110.32	159,00 €	7-14 days	622110.50	161,00 €	7-10 days
2,00 mm	Unperforated	185	3,30	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512108.32	151,00 €	7-14 days	512108.50	153,00 €	4-10 days
2,00 mm	Perforated	200	3,50	1,0 mm	Rv	ISO 3310-2 / ASTM E 323	542108.32	125,00 €	7-14 days	542108.50	127,00 €	4-10 days
2,24 mm	Unperforated	185	4,10	0,8 mm	Rv	ISO 3310-2	512109.32	156,00 €	7-14 days	512109.50	158,00 €	7-10 days
2,24 mm	Perforated	200	3,50	1,0 mm	Rv	ISO 3310-2	542110.32	125,00 €	7-14 days	542110.50	127,00 €	7-10 days
2,50 mm	Unperforated	185	4,00	0,8 mm	Rv	ISO 3310-2	512111.32	151,00 €	7-14 days	512111.50	153,00 €	7-10 days
2,50 mm	Perforated	200	4,00	1,0 mm	Rv	ISO 3310-2	542112.32	125,00 €	7-14 days	542112.50	127,00 €	7-10 days
2,58 mm	Unperforated	185	3,97	0,8 mm	Rv	U.S. CFR	622111.32	159,00 €	7-14 days	622111.50	161,00 €	7-10 days
2,80 mm	Unperforated	185	4,35	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512112.32	151,00 €	7-14 days	512112.50	153,00 €	7-10 days
3,00 mm	Unperforated	185	5,00	0,8 mm	Rv	Tolerances ISO 3310	520054.32	151,00 €	7-14 days	520054.50	153,00 €	7-10 days
3,00 mm	Perforated	200	5,00	1,0 mm	Rv	Tolerances ISO 3310	542116.32	125,00 €	7-14 days	542116.50	127,00 €	7-10 days
3,15 mm	Unperforated	185	4,70	0,8 mm	Rv	ISO 3310-2	512113.32	151,00 €	7-14 days	512113.50	153,00 €	4-10 days
3,25 mm	Unperforated	185	5,70	0,8 mm	Rv	Tolerances ISO 3310	520156.32	151,00 €	7-14 days	520156.50	153,00 €	7-10 days
3,55 mm	Unperforated	185	5,20	0,8 mm	Rv	ISO 3310-2	512115.32	151,00 €	7-14 days	512115.50	153,00 €	7-10 days
4,00 mm	Unperforated	185	5,80	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512116.32	151,00 €	7-14 days	512116.50	153,00 €	4-10 days
4,00 mm	Perforated	200	6,00	1,0 mm	Rv	Tolerances ISO 3310	542119.32	118,00 €	7-14 days	542119.50	120,00 €	4-10 days
4,25 mm	Unperforated	185	6,10	0,8 mm	Rv	Tolerances ISO 3310	520158.32	156,00 €	7-14 days	520158.50	158,00 €	7-10 days
4,50 mm	Unperforated	185	6,30	0,8 mm	Rv	ISO 3310-2	512117.32	151,00 €	7-14 days	512117.50	153,00 €	4-10 days
4,50 mm	Unperforated	185	6,30	0,8 mm	Rv	EN ISO 5223	732252.32	166,10 €	7-14 days	732252.50	168,10 €	7-10 days
4,75 mm	Unperforated	185	6,60	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512118.32	156,00 €	7-14 days	512118.50	158,00 €	4-10 days

**Perforated Plate Test Sieves with Round Apertures, Ø 200 mm, Stainless Steel** **Blau-Metall Specs & Prices (01/2025)**

Aperture (w)	Margin / Edge	Perf. Area (ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	200 x 32 mm			200 x 50 mm		
							Art-No.	€/Pcs.	Deliv. Time*	Art-No.	€/Pcs.	Deliv. Time*
4,76 mm	Unperforated	185	6,35	0,8 mm	Rv	U.S. CFR	622112.32	159,00 €	7-14 days	622112.50	161,00 €	7-10 days
4,80 mm	Unperforated	185	6,30	0,8 mm	Rv	Tolerances ISO 3310	520159.32	156,00 €	7-14 days	520159.50	158,00 €	7-10 days
5,00 mm	Unperforated	185	6,90	0,8 mm	Rv	ISO 3310-2	512119.32	151,00 €	7-14 days	512119.50	153,00 €	4-10 days
5,30 mm	Unperforated	185	8,10	1,0 mm	Rv	Tolerances ISO 3310	520160.32	174,00 €	7-14 days	520160.50	176,00 €	7-10 days
5,60 mm	Unperforated	185	7,70	0,8 mm	Rv	ISO 3310-2 / ASTM E 323	512120.32	151,00 €	7-14 days	512120.50	153,00 €	4-10 days
6,00 mm	Unperforated	185	8,00	0,8 mm	Rv	Tolerances ISO 3310	520056.32	151,00 €	7-14 days	520056.50	153,00 €	4-10 days
6,00 mm	Perforated	200	8,00	1,0 mm	Rv	Tolerances ISO 3310	542125.32	118,00 €	7-14 days	542125.50	120,00 €	7-10 days
6,30 mm	Unperforated	185	8,50	1,0 mm	Rv	ISO 3310-2 / ASTM E 323	512121.32	151,00 €	7-14 days	512121.50	153,00 €	7-10 days
6,50 mm	Unperforated	185	8,80	1,0 mm	Rv	Tolerances ISO 3310	520162.32	168,00 €	7-14 days	520162.50	170,00 €	7-10 days
6,70 mm	Unperforated	185	9,40	1,0 mm	Rv	ISO 3310-2	512122.32	151,00 €	7-14 days	512122.50	153,00 €	7-10 days
7,00 mm	Unperforated	185	9,40	0,8 mm	Rv	Tolerances ISO 3310	520057.32	151,00 €	7-14 days	520057.50	153,00 €	7-10 days
7,10 mm	Unperforated	185	9,40	1,0 mm	Rv	ISO 3310-2	512123.32	151,00 €	7-14 days	512123.50	153,00 €	4-10 days
7,50 mm	Unperforated	185	10,30	1,0 mm	Rv	Tolerances ISO 3310	520163.32	168,00 €	7-14 days	520163.50	170,00 €	7-10 days
8,00 mm	Unperforated	185	10,4	1,5 mm	Rv	ISO 3310-2 / ASTM E 323	512124.32	151,00 €	7-14 days	512124.50	153,00 €	4-10 days
8,00 mm	Perforated	200	10,0	1,0 mm	Rv	Tolerances ISO 3310	542132.32	118,00 €	7-14 days	542132.50	120,00 €	7-10 days
9,00 mm	Unperforated	185	11,0	1,5 mm	Rv	ISO 3310-2	512125.32	151,00 €	7-14 days	512125.50	153,00 €	7-10 days
9,50 mm	Unperforated	185	12,1	1,5 mm	Rv	ISO 3310-2 / ASTM E 323	512126.32	151,00 €	7-14 days	512126.50	153,00 €	4-10 days
10,0 mm	Unperforated	185	12,6	1,5 mm	Rv	ISO 3310-2	512127.32	151,00 €	7-14 days	512127.50	153,00 €	4-10 days
10,0 mm	Perforated	200	14,0	1,5 mm	Rv	Tolerances ISO 3310	542137.32	118,00 €	7-14 days	542137.50	120,00 €	7-10 days
11,0 mm	Unperforated	max. 185	15,0	1,5 mm	Rv	Tolerances ISO 3310	520166.32	168,00 €	7-14 days	520166.50	170,00 €	7-10 days
11,2 mm	Unperforated	max. 185	15,0	1,5 mm	Rv	ISO 3310-2 / ASTM E 323	512128.32	151,00 €	7-14 days	512128.50	153,00 €	4-10 days
12,0 mm	Unperforated	max. 185	15,5	1,5 mm	Rv	Tolerances ISO 3310	520168.32	168,00 €	7-14 days	520168.50	170,00 €	7-10 days
12,5 mm	Unperforated	max. 185	17,5	1,5 mm	Rv	ISO 3310-2 / ASTM E 323	512129.32	151,00 €	7-14 days	512129.50	153,00 €	7-10 days
13,0 mm	Unperforated	max. 185	17,5	1,5 mm	Rv	Tolerances ISO 3310	520170.32	168,00 €	7-14 days	520170.50	170,00 €	7-10 days
13,2 mm	Unperforated	max. 185	18,5	1,5 mm	Rv	ISO 3310-2 / ASTM E 323	512130.32	151,00 €	7-14 days	512130.50	153,00 €	7-10 days
14,0 mm	Unperforated	max. 185	18,5	1,5 mm	Rv	ISO 3310-2	512131.32	151,00 €	7-14 days	512131.50	153,00 €	7-10 days
15,0 mm	Unperforated	max. 185	20,0	2,0 mm	Rv	Tolerances ISO 3310	520172.32	168,00 €	7-14 days	520172.50	170,00 €	7-10 days
16,0 mm	Unperforated	max. 185	21,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512132.32	151,00 €	7-14 days	512132.50	153,00 €	4-10 days
17,0 mm	Unperforated	max. 185	21,5	2,0 mm	Rv	Tolerances ISO 3310	520173.32	168,00 €	7-14 days	520173.50	170,00 €	7-10 days
18,0 mm	Unperforated	max. 185	23,0	2,0 mm	Rv	ISO 3310-2	512133.32	151,00 €	14-21 days	512133.50	153,00 €	14-21 days
19,0 mm	Unperforated	max. 185	24,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512134.32	151,00 €	14-21 days	512134.50	153,00 €	14-21 days
20,0 mm	Unperforated	max. 185	25,0	2,0 mm	Rv	ISO 3310-2	512135.32	151,00 €	14-21 days	512135.50	153,00 €	14-21 days

# Perforated Plate Test Sieves with Round Apertures, Ø 200 mm, Stainless Steel

Blau-Metall Specs & Prices (01/2025)

Aperture (w)	Margin / Edge	Perf. Area (ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	200 x 32 mm			200 x 50 mm		
							Art-No.	€/Pcs.	Deliv. Time*	Art-No.	€/Pcs.	Deliv. Time*
22,4 mm	Unperforated	max. 185	28,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512136.32	151,00 €	7-14 days	512136.50	153,00 €	4-10 days
25,0 mm	Unperforated	max. 185	31,5	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512137.32	151,00 €	7-14 days	512137.50	153,00 €	7-14 days
26,5 mm	Unperforated	max. 185	34,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512138.32	151,00 €	14-21 days	512138.50	153,00 €	14-21 days
28,0 mm	Unperforated	max. 185	36,0	2,0 mm	Rv	ISO 3310-2	512139.32	151,00 €	7-14 days	512139.50	153,00 €	7-10 days
31,5 mm	Unperforated	max. 185	40,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512140.32	151,00 €	14-21 days	512140.50	153,00 €	14-21 days
35,5 mm	Unperforated	max. 185	45,0	2,0 mm	Rv	ISO 3310-2	512141.32	151,00 €	14-21 days	512141.50	153,00 €	14-21 days
37,5 mm	Unperforated	max. 185	47,5	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512142.32	151,00 €	14-21 days	512142.50	153,00 €	14-21 days
40,0 mm	Unperforated	max. 185	50,0	2,0 mm	Rv	ISO 3310-2	512143.32	151,00 €	14-21 days	512143.50	153,00 €	14-21 days
45,0 mm	Unperforated	max. 185	56,0	2,0 mm	Rv	ISO 3310-2 / ASTM E 323	512144.32	151,00 €	14-21 days	512144.50	153,00 €	14-21 days
50,0 mm	Unperforated	max. 185	63,0	3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512145.32	151,00 €	14-21 days	512145.50	153,00 €	14-21 days
53,0 mm	Unperforated	max. 185	67,0	3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512146.32	151,00 €	14-21 days	512146.50	153,00 €	14-21 days
56,0 mm	Unperforated	max. 185	71,0	3,0 mm	Rv	ISO 3310-2	512147.32	151,00 €	14-21 days	512147.50	153,00 €	14-21 days
63,0 mm	Unperforated	max. 185	80,0	3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512148.32	151,00 €	14-21 days	512148.50	153,00 €	14-21 days
71,0 mm	Unperforated	max. 185	90,0	3,0 mm	Rv	ISO 3310-2	512149.32	151,00 €	14-21 days	512149.50	153,00 €	14-21 days
75,0 mm	Unperforated	max. 185	95,0	3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512150.32	151,00 €	14-21 days	512150.50	153,00 €	14-21 days
80,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2	512151.32	151,00 €	14-21 days	512151.50	153,00 €	14-21 days
90,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512152.32	151,00 €	14-21 days	512152.50	153,00 €	14-21 days
100,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512153.32	151,00 €	14-21 days	512153.50	153,00 €	14-21 days
106,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512154.32	151,00 €	14-21 days	512154.50	153,00 €	14-21 days
112,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2	512155.32	151,00 €	14-21 days	512155.50	153,00 €	14-21 days
125,0 mm	Unperforated	Single Aperture		3,0 mm	Rv	ISO 3310-2 / ASTM E 323	512156.32	151,00 €	14-21 days	512156.50	153,00 €	14-21 days

**Other aperture sizes and/or tolerance specifications available on request.**

All sieves include a Certificate of Compliance 2.1 (in accordance with EN 10204) and a test sieve record card. (\*\* without documentation)

\* Delivery times are approximate and for domestic delivery within Germany. International delivery times vary country-to-country; express shipping options available.

\* Delivery times for bulk order of 15+ test sieves available on request.

## Optional Features

3D Laser Label and Sieve Documentation in English Language	--	0,00 €	7-14 days	--	0,00 €	7-14 days
Test Sieve according to Food Safety Standard (Solder joints of 99.9% Sn)	810230	+4,50 €	7-14 days	810230	+4,50 €	7-14 days
Temperature Resistance (up to 200°C) + Food-Safety	810231	+12,00 €	7-14 days	810231	+12,00 €	7-14 days
Inspection Certificate 3.1 (EN 10204, Certification acc. to ISO / ASTM)	910101	+22,50 €	--	910101	+22,50 €	--
Calibration Certificate 3.1 (EN 10204, Calibration acc. to ISO / ASTM)	910201	+40,00 €	--	910201	+40,00 €	--