

# Perforated Plate Test Sieves with Square Apertures, Ø 200 mm, Stainless Steel

# Blau-Metall Specs & Prices (10/2023)

Aperture (w)	U.S. Alternative	Perf. Area (Ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	200 x 32 mm			200 x 50 mm		
							Art-No.	€/Pcs.	Deliv. Time*	Art-No.	€/Pcs.	Deliv. Time*
4,00 mm	5/32 in.	max. 185	5,80	1,0 mm	Qg	ISO 3310-2 / ASTM E 323	513116.32	123,00 €	7-14 days	513116.50	125,00 €	3-7 days
4,50 mm		max. 185	6,30	1,0 mm	Qg	ISO 3310-2	513117.32	133,00 €	7-14 days	513117.50	135,00 €	4-7 days
4,75 mm	3/16 in.	max. 185	7,20	1,0 mm	Qg	ISO 3310-2 / ASTM E 323	513118.32	133,00 €	7-14 days	513118.50	135,00 €	4-7 days
5,00 mm		max. 185	6,90	1,0 mm	Qg	ISO 3310-2	513119.32	133,00 €	7-14 days	513119.50	135,00 €	4-7 days
5,60 mm	7/32 in.	max. 185	8,20/8,60	1,0 mm	Qg	ISO 3310-2 / ASTM E 323	513120.32	123,00 €	7-14 days	513120.50	125,00 €	3-7 days
6,00 mm		max. 185	9,00	1,0 mm	Qg	Tolerances (ISO 3310)	513210.32	138,00 €	7-14 days	513210.50	140,00 €	7-14 days
6,30 mm	¼ in.	max. 185	8,50	1,0 mm	Qg	ISO 3310-2 / ASTM E 323	513121.32	133,00 €	7-14 days	513121.50	135,00 €	4-7 days
6,70 mm	17/64 in.	max. 185	8,90	1,0 mm	Qg	ISO 3310-2 / ASTM E 323	513122.32	133,00 €	14-21 days	513122.50	135,00 €	14-21 days
7,00 mm		max. 185	10,2	1,5 mm	Qg	Tolerances (ISO 3310)	513214.32	138,00 €	14-21 days	513214.50	140,00 €	14-21 days
7,10 mm		max. 185	10,0	1,0 mm	Qg	ISO 3310-2	513123.32	133,00 €	7-14 days	513123.50	135,00 €	4-7 days
7,50 mm		max. 185	10,5	1,5 mm	Qg	Tolerances (ISO 3310)	513216.32	138,00 €	14-21 days	513216.50	140,00 €	14-21 days
8,00 mm	5/16 in.	max. 185	10,4	1,5 mm	Qg	ISO 3310-2 / ASTM E 323	513124.32	123,00 €	7-14 days	513124.50	125,00 €	3-7 days
9,00 mm		max. 185	12,0	1,5 mm	Qg	ISO 3310-2	513125.32	133,00 €	7-14 days	513125.50	135,00 €	4-7 days
9,50 mm	¾ in.	max. 185	13,1	1,5 mm	Qg	ISO 3310-2 / ASTM E 323	513126.32	133,00 €	7-14 days	513126.50	135,00 €	4-7 days
10,0 mm		max. 185	12,6	1,5 mm	Qg	ISO 3310-2	513127.32	123,00 €	7-14 days	513127.50	125,00 €	4-7 days
11,0 mm		max. 185	15,0	1,5 mm	Qg	Tolerances (ISO 3310)	513224.32	138,00 €	14-21 days	513224.50	140,00 €	14-21 days
11,2 mm	7/16 in.	max. 185	15,0	1,5 mm	Qg	ISO 3310-2 / ASTM E 323	513128.32	123,00 €	7-14 days	513128.50	125,00 €	3-7 days
12,0 mm		max. 185	16,0	1,5 mm	Qg	Tolerances (ISO 3310)	513227.32	138,00 €	14-21 days	513227.50	140,00 €	14-21 days
12,5 mm	½ in.	max. 185	17,5	1,5 mm	Qg	ISO 3310-2 / ASTM E 323	513129.32	133,00 €	7-14 days	513129.50	135,00 €	4-7 days
13,0 mm		max. 185	18,0	1,5 mm	Qg	Tolerances (ISO 3310)	513230.32	138,00 €	14-21 days	513230.50	140,00 €	14-21 days
13,2 mm	17/32 in.	max. 185	18,5	1,5 mm	Qg	ISO 3310-2 / ASTM E 323	513130.32	133,00 €	7-14 days	513130.50	135,00 €	4-7 days
14,0 mm		max. 185	18,0	1,5 mm	Qg	ISO 3310-2	513131.32	133,00 €	7-14 days	513131.50	135,00 €	7-14 Tage
15,0 mm		max. 185	19,0	2,0 mm	Qg	Tolerances (ISO 3310)	513233.32	138,00 €	14-21 days	513233.50	140,00 €	14-21 days
16,0 mm	5/8 in.	max. 185	20,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513132.32	123,00 €	7-14 days	513132.50	125,00 €	3-7 days
18,0 mm		max. 185	23,0	2,0 mm	Qg	ISO 3310-2	513133.32	133,00 €	7-14 days	513133.50	135,00 €	7-14 days
19,0 mm	¾ in.	max. 185	24,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513134.32	133,00 €	7-14 days	513134.50	135,00 €	7-14 days
20,0 mm		max. 185	25,0	2,0 mm	Qg	ISO 3310-2	513135.32	133,00 €	7-14 days	513135.50	135,00 €	7-14 days
22,0 mm		max. 185	27,0	2,0 mm	Qg	Tolerances (ISO 3310)	513244.32	138,00 €	14-21 days	513244.50	140,00 €	14-21 days
22,4 mm	⅞ in.	max. 185	28,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513136.32	123,00 €	7-14 days	513136.50	125,00 €	3-7 days
24,0 mm		max. 185	29,0	2,0 mm	Qg	Tolerances (ISO 3310)	513248.32	138,00 €	14-21 days	513248.50	140,00 €	14-21 days
25,0 mm	1 in.	max. 185	31,5	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513137.32	128,00 €	7-14 days	513137.50	130,00 €	7-14 Tage
26,0 mm		max. 185	32,0	2,0 mm	Qg	Tolerances (ISO 3310)	513251.32	138,00 €	14-21 days	513251.50	140,00 €	14-21 days
26,5 mm	1 1/16 in.	max. 185	34,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513138.32	128,00 €	7-14 days	513138.50	130,00 €	7-14 Tage

# Perforated Plate Test Sieves with Square Apertures, Ø 200 mm, Stainless Steel

# Blau-Metall Specs & Prices (10/2023)

Aperture (w)	U.S. Alternative	Perf. Area (Ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	200 x 32 mm			200 x 50 mm		
							Art-No.	€/Pcs.	Deliv. Time*	Art-No.	€/Pcs.	Deliv. Time*
28,0 mm		max. 185	36,0	2,0 mm	Qg	ISO 3310-2	513139.32	128,00 €	7-14 days	513139.50	130,00 €	7-14 Tage
30,0 mm		max. 185	37,0	2,0 mm	Qg	Tolerances (ISO 3310)	513257.32	138,00 €	14-21 days	513257.50	140,00 €	14-21 days
31,5 mm	1¼ in.	max. 185	40,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513140.32	113,00 €	7-14 days	513140.50	115,00 €	3-7 days
35,5 mm		max. 185	45,0	2,0 mm	Qg	ISO 3310-2	513141.32	128,00 €	7-14 days	513141.50	130,00 €	7-14 days
37,5 mm	1½ in.	max. 185	47,5	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513142.32	128,00 €	7-14 days	513142.50	130,00 €	7-14 days
40,0 mm		max. 185	50,0	2,0 mm	Qg	ISO 3310-2	513143.32	128,00 €	7-14 days	513143.50	130,00 €	7-14 days
41,5 mm		max. 185	51,0	2,0 mm	Qg	Tolerances (ISO 3310)	513263.32	138,00 €	7-14 days	513263.50	140,00 €	7-14 days
45,0 mm	1¾ in.	max. 185	56,0	2,0 mm	Qg	ISO 3310-2 / ASTM E 323	513144.32	113,00 €	7-14 days	513144.50	115,00 €	3-7 days
50,0 mm	2 in.	max. 185	63,0	3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513145.32	128,00 €	14-21 days	513145.50	130,00 €	14-21 days
53,0 mm	2½ in.	max. 185	67,0	3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513146.32	128,00 €	14-21 days	513146.50	130,00 €	14-21 days
56,0 mm		max. 185	71,0	3,0 mm	Qg	ISO 3310-2	513147.32	128,00 €	7-14 days	513147.50	130,00 €	4-7 days
63,0 mm	2½ in.	max. 185	80,0	3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513148.32	113,00 €	7-14 days	513148.50	115,00 €	3-7 days
71,0 mm		Single Aperture		3,0 mm	Qg	ISO 3310-2	513149.32	113,00 €	14-21 days	513149.50	115,00 €	14-21 days
75,0 mm	3 in.	Single Aperture		3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513150.32	113,00 €	14-21 days	513150.50	115,00 €	14-21 days
80,0 mm		Single Aperture		3,0 mm	Qg	ISO 3310-2	513151.32	113,00 €	14-21 days	513151.50	115,00 €	14-21 days
90,0 mm	3½ in.	Single Aperture		3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513152.32	113,00 €	14-21 days	513152.50	115,00 €	14-21 days
100,0 mm	4 in.	Single Aperture		3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513153.32	113,00 €	14-21 days	513153.50	115,00 €	14-21 days
106,0 mm	4¼ in.	Single Aperture		3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513154.32	113,00 €	14-21 days	513154.50	115,00 €	14-21 days
112,0 mm		Single Aperture		3,0 mm	Qg	ISO 3310-2	513155.32	113,00 €	14-21 days	513155.50	115,00 €	14-21 days
125,0 mm	5 in.	Single Aperture		3,0 mm	Qg	ISO 3310-2 / ASTM E 323	513156.32	113,00 €	14-21 days	513156.50	115,00 €	14-21 days
Custom Aperture up to 10,0 mm				1,0-1,5 mm	Qg/Qv	Tolerance values of ISO 3310	--	176,00 €	14-21 days	--	178,00 €	14-21 days
Custom Aperture up to 22,4 mm				1,5-2,0 mm	Qg/Qv	Tolerance values of ISO 3310	--	158,00 €	14-21 days	--	160,00 €	14-21 days
Custom Aperture up to 125 mm				1,5-3,0 mm	Qg/Qv	Tolerance values of ISO 3310	--	144,00 €	14-21 days	--	146,00 €	14-21 days

All sieves include a Certificate of Compliance 2.1 (in accordance with EN 10204) and a test sieve record card.

\* Delivery times are approximate and for domestic delivery within Germany. International delivery times vary country-to-country; express shipping options available.

\* Delivery times for bulk order of 15+ test sieves available on request.

## Optional Features

3D Laser Label and Sieve Documentation in English Language	--	0,00 €	10-21 days	--	0,00 €	10-21 days
Test Sieve according to Food Safety Standard (Solder joints of 99.9% Sn)	810230	+4,50 €	10-21 days	810230	+4,50 €	10-21 days
Temperature Resistance (up to 200°C) + Food-Safety	810231	+12,00 €	10-21 days	810231	+12,00 €	10-21 days
Inspection Certificate 3.1 (EN 10204, Certification acc. to ISO / ASTM)	910101	+20,00 €	--	910101	+20,00 €	--
Calibration Certificate 3.1 (EN 10204, Calibration acc. to ISO / ASTM)	910201	+38,00 €	--	910201	+38,00 €	--