

Perforated Plate Sieves with Oblong Apertures, Ø 300 mm, Stainless Steel **Blau-Metall Specs & Prices (10/2023)**

Aperture (w ₁ x w ₂)	Tolerance (Δw ₁ x Δw ₂)	Pitch (p ₁ x p ₂)	Tolerance (Δp ₁ x Δp ₂)	Perf. Area (ø/mm)	Plate Thickness	Perf. Type	Standard	300 x 65 mm		
								Art-No.	€/Pcs.	Deliv. Time*
1,00 x 20,0	±0,03 x ±0,20	3,00 x 25,0	±0,20 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832001	265,00 €	~ 4 w.
1,10 x 20,0	+0,03/-0,06 x ±0,20	2,50 x 24,0	±0,20 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832110	305,00 €	~ 4 w.
1,20 x 20,0	±0,03 x ±0,20	2,60 x 24,0	±0,20 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832111	305,00 €	~ 4 w.
1,30 x 20,0	±0,03 x ±0,20	2,80 x 24,0	±0,20 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832112	305,00 €	~ 4 w.
1,40 x 20,0	+0,03/-0,05 x ±0,20	3,00 x 24,0	±0,20 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832113	305,00 €	~ 4 w.
1,50 x 20,0	±0,04 x ±0,20	4,00 x 25,0	±0,24 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832002	265,00 €	~ 4 w.
1,60 x 20,0	±0,04 x ±0,20	4,00 x 25,0	±0,24 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832003	265,00 €	~ 4 w.
1,70 x 20,0	±0,04 x ±0,20	4,00 x 25,0	±0,24 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832004	265,00 €	~ 4 w.
1,80 x 20,0	±0,04 x ±0,20	4,20 x 25,0	±0,24 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832005	265,00 €	~ 4 w.
1,90 x 20,0	±0,04 x ±0,20	4,30 x 25,0	±0,24 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832006	265,00 €	~ 4 w.
2,00 x 20,0	±0,04 x ±0,20	4,50 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832007	265,00 €	~ 4 w.
2,20 x 20,0	±0,05 x ±0,20	4,90 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832008	265,00 €	~ 4 w.
2,25 x 20,0	±0,05 x ±0,20	4,90 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832009	265,00 €	~ 4 w.
2,50 x 20,0	±0,05 x ±0,20	4,90 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832010	265,00 €	~ 4 w.
2,80 x 20,0	±0,05 x ±0,20	4,90 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832011	265,00 €	~ 4 w.
3,00 x 20,0	±0,05 x ±0,20	6,30 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832012	285,00 €	~ 4 w.
3,25 x 20,0	±0,05 x ±0,20	6,50 x 25,0	±0,26 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832013	285,00 €	~ 4 w.
3,50 x 20,0	±0,06 x ±0,20	6,80 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832014	265,00 €	~ 4 w.
3,55 x 20,0	±0,06 x ±0,20	6,80 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832015	265,00 €	~ 4 w.
3,75 x 20,0	±0,07 x ±0,20	7,25 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832116	285,00 €	~ 4 w.
4,00 x 20,0	±0,07 x ±0,20	7,40 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832117	265,00 €	~ 4 w.
4,25 x 20,0	±0,07 x ±0,20	7,70 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832118	265,00 €	~ 4 w.
4,50 x 20,0	±0,07 x ±0,20	8,00 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832119	265,00 €	~ 4 w.
4,75 x 20,0	±0,07 x ±0,20	8,30 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832120	265,00 €	~ 4 w.
5,00 x 20,0	±0,08 x ±0,20	9,00 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832121	265,00 €	~ 4 w.
6,00 x 20,0	±0,08 x ±0,20	10,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832125	265,00 €	~ 4 w.
7,00 x 20,0	±0,08 x ±0,20	11,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832129	265,00 €	~ 4 w.
7,50 x 20,0	±0,08 x ±0,20	12,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832131	265,00 €	~ 4 w.
8,00 x 20,0	±0,08 x ±0,20	12,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832133	265,00 €	~ 4 w.
8,00 x 12,0	±0,08 x ±0,20	12,5 x 17,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832220	295,00 €	~ 4 w.
10,0 x 20,0	±0,08 x ±0,20	14,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832141	265,00 €	~ 4 w.
14,0 x 20,0	±0,08 x ±0,20	18,0 x 25,0	±0,34 x ±0,50	max. 285	1,0 mm	Lg	Tolerance ISO 5223	832145	265,00 €	~ 4 w.

Other aperture sizes and/or tolerance specifications available on request.

All sieves include a Certificate of Compliance 2.1 (in accordance with EN 10204) and a test sieve record card.

* Delivery times are approximate and for domestic delivery within Germany. International delivery times vary country-to-country; express shipping options available. Delivery times for bulk order of 15+ test sieves available on request.

Optional Features

3D Laser Label and Sieve Documentation in English Language	--	0,00 €	~ 4 w.
2x Stainless Steel Handles for Sieving by Hand (Handles welded to sieve frame)	811205	+35,00 €	~ 5 w.
Test Sieve according to Food Safety Standard (Solder joints of 99.9% Sn)	811230	+21,00 €	~ 4 w.
Temperature Resistance (up to 200°C) + Food-Safety	811231	+40,00 €	~ 4 w.
Inspection Certificate 3.1 (EN 10204, Certification acc. to ISO / ASTM)	911101	+35,00 €	--
Calibration Certificate 3.1 (EN 10204, Calibration acc. to ISO / ASTM)	911201	+55,00 €	--