

Aperture (w)	U.S. Alternative	Perf. Area (Ø/mm)	Pitch (p/mm)	Plate Thickness	Perf. Type	Standard	300 x 65 mm		
							Art-No.	€/Pcs.	Deliv. Time*
4,00 mm	5/32 in.	max. 285	5,80	1,0 mm	Qg	ISO 3310-2/ASTM E 323	613116	163,00 €	7-14 days
4,50 mm		max. 285	7,00	1,0 mm	Qg	ISO 3310-2	613117	183,00 €	On Request
4,75 mm	3/16 in.	max. 285	7,30	1,0 mm	Qg	ISO 3310-2/ASTM E 323	613118	183,00 €	On Request
5,00 mm		max. 285	7,90	1,0 mm	Qg	ISO 3310-2	613119	183,00 €	On Request
5,60 mm	7/32 in.	max. 285	8,60	1,0 mm	Qg	ISO 3310-2/ASTM E 323	613120	178,00 €	14-21 days
6,30 mm	¼ in.	max. 285	9,00	1,0 mm	Qg	ISO 3310-2/ASTM E 323	613121	183,00 €	14-21 days
6,70 mm	17/64 in.	max. 285	9,00	1,0 mm	Qg	ISO 3310-2/ASTM E 323	613122	183,00 €	14-21 days
7,10 mm		max. 285	10,0	1,5 mm	Qg	ISO 3310-2	613123	183,00 €	14-21 days
8,00 mm	5/16 in.	max. 285	10,4	1,5 mm	Qg	ISO 3310-2/ASTM E 323	613124	163,00 €	7-14 days
9,00 mm		max. 285	13,0	1,5 mm	Qg	ISO 3310-2	613125	183,00 €	14-21 days
9,50 mm	⅜ in.	max. 285	13,2	1,5 mm	Qg	ISO 3310-2/ASTM E 323	613126	183,00 €	14-21 days
10,0 mm		max. 285	13,5	1,5 mm	Qg	ISO 3310-2	613127	183,00 €	14-21 days
11,2 mm	7/16 in.	max. 285	14,0	1,5 mm	Qg	ISO 3310-2/ASTM E 323	613128	163,00 €	7-14 days
12,5 mm	½ in.	max. 285	16,0	1,5 mm	Qg	ISO 3310-2/ASTM E 323	613129	183,00 €	14-21 days
13,2 mm	17/32 in.	max. 285	17,0	1,5 mm	Qg	ISO 3310-2/ASTM E 323	613130	183,00 €	14-21 days
14,0 mm		max. 285	18,0	1,5 mm	Qg	ISO 3310-2	613131	183,00 €	14-21 days
16,0 mm	5/8 in.	max. 285	20,0	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613132	163,00 €	7-14 days
18,0 mm		max. 285	22,4	2,0 mm	Qg	ISO 3310-2	613133	178,00 €	14-21 days
19,0 mm	¾ in.	max. 285	23,6	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613134	178,00 €	14-21 days
20,0 mm		max. 285	25,0	2,0 mm	Qg	ISO 3310-2	613135	178,00 €	14-21 days
22,4 mm	⅞ in.	max. 285	28,0	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613136	163,00 €	7-14 days
25,0 mm	1 in.	max. 285	31,5	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613137	178,00 €	14-21 days
26,5 mm	1 1/16 in.	max. 285	33,5	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613138	178,00 €	14-21 days
28,0 mm		max. 285	35,5	2,0 mm	Qg	ISO 3310-2	613139	178,00 €	14-21 days
31,5 mm	1¼ in.	max. 285	40,0	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613140	163,00 €	7-14 days
35,5 mm		max. 285	45,0	2,0 mm	Qg	ISO 3310-2	613141	178,00 €	14-21 days
37,5 mm	1½ in.	max. 285	47,5	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613142	178,00 €	14-21 days
40,0 mm		max. 285	50,0	2,0 mm	Qg	ISO 3310-2	613143	178,00 €	14-21 days
45,0 mm	1¾ in.	max. 285	56,0	2,0 mm	Qg	ISO 3310-2/ASTM E 323	613144	163,00 €	7-14 days
50,0 mm	2 in.	max. 285	63,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613145	178,00 €	14-21 days
53,0 mm	2⅝ in.	max. 285	67,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613146	178,00 €	14-21 days
56,0 mm		max. 285	71,0	3,0 mm	Qg	ISO 3310-2	613147	178,00 €	7-14 days
63,0 mm	2½ in.	max. 285	82,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613148	158,00 €	7-14 days
71,0 mm		max. 285	90,0	3,0 mm	Qg	ISO 3310-2	613149	168,00 €	14-21 days
75,0 mm	3 in.	max. 285	95,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613150	168,00 €	14-21 days
80,0 mm		max. 285	100,0	3,0 mm	Qg	ISO 3310-2	613151	168,00 €	14-21 days
90,0 mm	3½ in.	max. 285	112,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613152	158,00 €	14-21 days
100,0 mm	4 in.	max. 285	125,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613153	168,00 €	14-21 days
106,0 mm	4¼ in.	max. 285	132,0	3,0 mm	Qg	ISO 3310-2/ASTM E 323	613154	168,00 €	14-21 days
112,0 mm		max. 285	141,0	3,0 mm	Qg	ISO 3310-2	613155	168,00 €	14-21 days
125,0 mm	5 in.	Single Aperture		3,0 mm	Qg	ISO 3310-2/ASTM E 323	613156	158,00 €	14-21 days
Custom Aperture up to 10,0 mm				1,0-1,5 mm	Qg/Qv	Tolerance values of ISO 3310	--	208,00 €	On Request
Custom Aperture up to 22,4 mm				1,5-2,0 mm	Qg/Qv	Tolerance values of ISO 3310	--	198,00 €	14-21 days
Custom Aperture up to 125 mm				1,5-3,0 mm	Qg/Qv	Tolerance values of ISO 3310	--	188,00 €	14-21 days

All sieves include a Certificate of Compliance 2.1 (in accordance with EN 10204) and a test sieve record card.

* Delivery times are approximate and for domestic delivery within Germany. International delivery times vary country-to-country; express shipping options available. Delivery times for bulk order of 15+ test sieves available on request.

Optional Features

3D Laser Label and Sieve Documentation in English Language	--	0,00 €	14-21 days
Test Sieve according to Food Safety Standard (Solder joints of 99.9% Sn)	811230	+21,00 €	14-21 days
Temperature Resistance (up to 200°C) + Food-Safety	811231	+30,00 €	14-21 days
Inspection Certificate 3.1 (EN 10204, Certification acc. to ISO / ASTM)	911101	+35,00 €	--
Calibration Certificate 3.1 (EN 10204, Calibration acc. to ISO / ASTM)	911201	+55,00 €	--